

SYSTEMATIC DECUBITUS THERAPY



POSITIONING AIDS **POSIMED®**



Funke
MEDICAL

Funke-Medical AG
Ährenfeld 10
46348 Raesfeld
Tel.: +49-2865-6648
Fax: +49-2865-8133
info@funke-medical.de
www.funke-medical.de

Positioning aids

Positioning aids in various forms and sizes for optimum positioning as a means of decubitus prophylaxis and treatment support.

With the various positioning aids the patient is moved marginally and the pressure distribution of the main support points enables the patient to lie in a relaxed and pressure-relieving manner.

The aim is to keep the duration and intensity of the effect of pressure as small as possible.

Our positioning system is based on a pressure-relieving and positionable filling of small spheres!

Advantages

- waterproof, but breathable
- virus and bacteria-proof
- disinfected by wiping
- autoclavable
- food-safe
- very low own weight
- rapid optimal adaptation to the contour of the body
- optimal microclimate
- optimal pressure relief → suitable for decubitus therapy
- 5-year guarantee on functional properties



Flame retardant according to DIN EN 597 (1+2)



Waterproof but breathable



CE marking



Virus and bacteria-proof



Do not wash



Autoclavable



Universal cushion

- For bedding and positioning

Size

37/26

56/40

Art. - Nr.

PO-UNIK3726

PO-UNIG5640



Heel rail

- For positioning the heel

Size
39/21/8

Art. - Nr.
PO-FERS3921



Forearm splint

- For the positioning of the forearm

Size
39/21/8

Art. - Nr.
PO-UNTE3921



Bolster

- For supporting and positioning the knee, neck and other extremities

Size
60/21

Art. - Nr.
PO-ROLL6021



Abduction pillow

- For the positioning of the knee, ankle and epicondyle


Size
60/35

Art. - Nr.
PO-SPRE6035



30° positioning cushion

- For providing support in the stable 30° position, with integrated pillow



Size
170/74
182/54

Art. - Nr.
PO-LAGE0165
PO-LAGE1165

Hip abduction pillow

- For the positioning of the hips, reduces the danger of dislocating the hip prostheses

Size
35/28/14

Art. - Nr.
PO-HÜFT3528



Semi-circle

- For positioning the upper body


Size
137/80

Art. - Nr.
PO-HALB1378




Positioning examples


30° lateral position



heel altitude position



half sitting



135° position

